


HydroCision™

TenJet™

**BRIDGING THE
GAP IN
TENDINOSIS
TREATMENT OPTIONS**



TenJet addresses the pathology, not just the symptoms

High velocity saline jet, with pressures up to 14,000 psi, acts as a cutting blade to resect degenerative tendon tissue without generating heat.


Venturi effect creates localized suction and pulls the degenerative, frayed, tissue inside a tiny, **1.65 mm**, cutting window while leaving healthy tissue unharmed

Instantaneous evacuation of diseased tissue allows for real-time ultrasound guided assessment of treatment


12-gauge 3-inch needle enables minimally invasive resection in all joints, including the challenging hip and shoulder anatomy



Minimally invasive treatment for tendinosis



Pre Operative




Post Operative 6 months

Achilles Tendinosis

60-year old patient had mid-substance and insertional Achilles tendinosis. He had been in chronic pain for 9 months with no improvement from PT. Post-TenJet, the patient was pain free at 3 months, with ultrasound showing complete remodeling at 6 months.




Pre Operative




Post Operative 3 months

Calcific Shoulder Tendinopathy

61-year old male presented with a painful right shoulder. Pre-op radiographs showed a calcific deposit in the supraspinatus. The patient had failed 9 months of conservative treatment. 12-week post-op radiographs revealed complete resolution of the deposit. Patient was pain free and had complete ROM.



Pre Operative



Post Operative 6 weeks

Medial Elbow Tendinosis

55-year old insulin dependent golfer with 2 years of medial elbow pain had been through PT program and received 1 CSI. Pre-procedure VAS score was 6-8 with activity. Post-TenJet, patient was pain-free within 3 weeks. 6-week ultrasound shows resolution of calcifications and reduction in tendinosis with contraction of the tendon.

Transform patient care with TenJet


To learn more, visit our product page at www.TenJet.com

Ordering Information

The TenJet device is used in conjunction with the HydroCision console.

TenJet kit includes a single use, sterile device and an accessory kit.

The accessory kit includes Chloraprep™, ultrasound probe cover, ultrasound gel, a 10 cc syringe, two 23-gauge needles, gauzes, Steri-strips™ and Tegaderm™.



Reference #	Description
52700	HydroCision Console
76000	TenJet, single unit

HydroCision Inc.

267 Boston Road, Suite 28
North Billerica, MA 01862

www.hydrocision.com

Telephone: 1-(978) 474 9300

Email: info@hydrocision.com

Orders/Inquiries: 1-(888) 747 4470

Chloraprep is a trademark of Carefusion 2200, Inc.

Steri-strips and Tegaderm are trademarks of 3M Company.

TenJet and HydroCision are trademarks of HydroCision Inc.

HydroCision and TenJet logo are trademarks of HydroCision Inc.

© 2020 HydroCision, Inc.

All rights reserved. 6000-0061C